
DE Digital Thermometer BigScreen #2142

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
EN Digital Thermometer BigScreen #2142

Please read the instructions carefully before using the thermometer for the first time, and keep these in a safe place.

Correct use

The fever thermometer is designed exclusively for use in the veterinary sector and not for use on humans.

Use

Before using, please disinfect the probe at first. To switch on, press the ON/OFF button next to the display; a short beep will sound, indicating that the thermometer is operational. When the letters „Lo“ or „Lo°F“ will appear on the LCD and if it is more than 42,7°C or 109,2°F, then „Hi°C“ or „Hi°F“ will appear on the LCD. Temperature below 37,6°C = Normal alarm: green backlight, 5 sec. Temperature above 37,6°C = Fever alarm: red backlight, 5 sec. How to change „Lo“ to „Lo°F“. Press switch to turn on and hold it, until displays „Lo°C“ (no flash), then release the switch within 2 sec. and press it again. You can change the temperature unit from °C to °F or back. As soon as the constant temperature value is reached, a beep will sound four times, and the „°C“ or „°F“ symbol will stop flashing. The highest temperature measured appears on the LCD. When the measurement is completed, please switch the thermometer off by pressing the ON/OFF button. After the temperature has been displayed, the thermometer will shut off automatically after 10 minutes.

Memory function

Switch the thermometer off, then press the ON/OFF button for about 2 seconds. The last measured value with „°C“ or „°F“ will appear automatically on the LCD.

Methods of measuring temperature In the rectum (rectal)

This is the most accurate method from a medical point of view. The thermometer tip is inserted carefully into the rectum for a maximum of 2 cm. The usual measuring time is approximately 9 to 11 seconds.

Change the measuring unit

To switch the unit between °C and °F, turn the thermometer off. Press and hold the „ON/OFF“ button until the LCD shows as Fig. A, then release the button. And press the button to choose the unit, the chosen unit will flashing. The thermometer will shut down automatically in 5 seconds.

Cleaning and disinfection

The best way to clean the thermometer tip is by applying a disinfectant (e.g. 70% ethyl alcohol) with a damp cloth. This thermometer is IP22 water resistant and can be protected against vertically falling water drops when enclosure tilted up to 15 degree through cleaning and disinfection.

Safety precautions

- Do not allow the device to come into contact with hot water.
- Do not expose to high temperatures or direct sunlight.
- Do not drop the thermometer. It is neither shock-proof nor impact-resistant.
- Do not bend or open the device (except the battery compartment).
- Do not clean with thinners, petrol or benzene. Only clean with water or disinfectant.
- Do not immerse the Waterproof thermometers in liquid for longer than 30 minutes.

- The thermometer contains small parts (battery, battery compartment) which can be swallowed by children. For this reason, do not leave the thermometer unattended in the hands of children.
- Avoid bending the Non-Flexible thermometer tip. For the Flexible thermometer, avoid bending the tip by more than 45 degrees.
- If the ambient temperature is over 35 °C or 95 °F, dip the thermometer tip in cold water for approx. 5 to 10 seconds prior to measuring the temperature.
- A persistent fever requires medical attention - please consult a veterinarian!

Battery replacement

If the battery icon  is shown on the LCD display, the thermometer is powered adequately. If the icon no longer appears or if the device will not switch on, the battery is empty and must be replaced. To do this, remove the battery compartment cover and replace the battery with a new one of the same type. Ensure the battery is inserted correctly. (+ top, - bottom).

Technical data

Measurement range: 32,0 ~ 42,9 °C / 89,6 ~ 109,2 °F
Measurement accuracy: ± 0,1 °C / 0,2 °F
(35,5 °C ~ 42,0 °C / 95,9 °F ~ 107,6 °F)
± 0,2 °C / 0,4 °F
(32,0 °C ~ 35,5 °C, 42,0 °C ~ 42,9 °C / 89,6 °F ~ 95,9 °F, 107,6 °F ~ 109,2 °F)

Storage/transportation

temperature:
Ambient temperature during use: (5~35) °C ≤ 80%RH
Min Scale: 0,1 °C / 0,1 °F
Atmospheric pressure: 700 ~ 1060 hPa
CR2032 Lithium, 3.0 V, lifetimium
100 hours under continuous operation
Weight: Approx. 27 g

Calibration check

This thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the operation instruction, periodic re-adjustment is not required.

CE-UKCA-conformity declaration


Albert Kerbl GmbH hereby declares that the product / device described in these instructions complies with the fundamental requirements and other relevant stipulations and regulations. The CE-UKCA-mark confirms compliance with the Directives of the European Union or the relevant UK legislation.

Electronic scrap

Do not dispose of batteries in domestic waste. Dispose of the battery properly, separately from the electrical device and in accordance with locally applicable regulations. Municipal collection points, as well as return points operated by retailers, can be used to return such batteries. The separate collection and reprocessing of old batteries is intended to avoid negative impacts on the environment and human health.

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Albert Kerbl GmbH Felizenzell 9
Kerbl UK Ltd Lands End Way Oakham, Rutland